

CLINEB is a professional piston compressor nebulizer for aerosoltherapy. Designed for intensive treatments in clinics or hospitals with large compartment for convenient storage of accessories and medications. Provided with pressure gauge and indicator.



Clineb ref.:

**RE 300550**

#### STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
  - Adult Masks with elastic strings
  - Child Masks with elastic strings
  - Mouthpiece
  - Nosepiece (non invasive)
  - Air Tube
  - Air Filter
- For spare consumables see page 27*



#### OPTIONS

- Stand on castors (ref. 27731)

#### TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	184 VA
Max Pressure	2.5 bar 36 psi 250 kPa
Operating Pressure	1.30 bar 18.9 psi 130 kPa
Max Air Flow	16 l/min
Operating Air Flow	5,2 l/min
Neb-rate with 4 ml solution	0,50 ml/min
Mass Median Aerodynamic Diameter	3,25 µm (MMAD) according to EN-13544-1
Respirable Fraction (< 5µm)	> 70%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 57 dBA
Weight	2,20 kg
Size (Main Unit)	cm 25 x 19 x 16
Years of Warranty	2
Place of Manufacturing	Italy
Shipping carton	4

**CE 0123**