

FLO-Baby is a compact and handy ultrasonic aerosol delivery system. It is recommended for infants, children and adults. Its absolute quietness and fastness in delivering medications make it a versatile portable nebulizer. Light-weight, easy to operate, with integrated mouthpiece and a support base for use on the table. It can be used as "hand-held" or "desk-type" unit depending on patient's convenience. Automatic shut-off system after ten minutes. In case of poor cooling or unintentional start of the device. Three different air flow settings to adjust treatment efficiency.



Flo Baby ref.:

RE 300900/02

STANDARD CONSUMABLES INCLUDED

- Main Unit with integrated mouthpiece and AC/DC Transformer
 - Adult Masks with elastic strings
 - Puppet-shaped Infant Mask with Fruit Fragrance
 - Corrugated Tube
 - Unit Base Support for Use on the Table
 - Air Filters (x2)
 - Medication Cups (x10)
 - Carrying Bag for Product Storage and Transport
- For spare consumables see page 27*

OPTIONS

- Rechargeable battery pack

CE 0123

TECHNICAL FEATURES

Transducer operating frequency	2,5MHz+200kHz (-100 kHz)
Power Feeding	230V-50Hz by AC/DC Transformer (included) 12V DC by Rechargeable Battery (Optional)
Power Consumption (AC/DC adapter)	10 VA
Noise Level	40 dBA
Medication cup max capacity	7.5 ml
Automatic power shut-off	11 min (+/- 2 min)
Nebulization rate	max 0.7ml/min with saline solution
Particle size	1-10µm
Size	cm 6 x 11x 16
Weight	g 220
Years of Warranty	2
Shipping Carton	6

