

EOLO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with convenient handle and a large compartment for easy transport of unit, accessories and medications. Recently upgraded for non-stop operation in case intensive treatments are required.



EOLO ref.:

**RE 300400**

#### STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
  - Adult Masks with elastic strings
  - Child Masks with elastic strings
  - Mouthpiece
  - Nosepiece (non invasive)
  - Air Tube
  - Air Filter
- For spare consumables see page 27*

#### TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	170 VA
Max Pressure	2.5 bar 36 psi 250 kPa
Operating Pressure	1.1 bar 16 psi 110 kPa
Max Air Flow	14 l/min
Operating Air Flow	5 l/min
Neb-rate with 4 ml solution	0,40 ml/min
Mass Median Aerodynamic Diameter	2,44 µm (MMAD) according to EN-13544-1
Respirable Fraction (< 5µm)	> 80%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 55 dBA
Weight	1,65 kg
Size (Main Unit)	cm 14 x 12 x 22
Years of Warranty	2
Place of Manufacturing	Italy
Shipping carton	6

**CE 0123**