

LITE is a miniaturized micropiston compressor system for aerosol drug delivery featuring highly efficient nebulizer. Noise level and vibrations are minimized thanks to the latest software "VRS" Vibration Reduction System. Small in size, easy to carry, compact and simple to use, LITE provides high quality of aerosol with any medication and it is therefore recommended for home-care and hospital use, for both adults and children whenever PORTABILITY is required.

Lite ref.:

RE 300850



STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
 - Adult Masks with elastic strings
 - Child Masks with elastic strings
 - Mouthpiece
 - Nosepiece (non invasive)
 - Air Tube
 - Air Filter
 - Carrying Bag
 - Universal Switching Adapter
- For spare consumables see page 27*



TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	100-240V – 50/60Hz
Power consumption	38,5 VA
Max Pressure	2 bar 29 psi 200 kPa
Operating Pressure	0,8 bar 11,6 psi 80 kPa
Max Air Flow	12 l/min
Operating Air Flow	4 l/min
Neb-rate with 4 ml solution	0,35 ml/min
Mass Median Aerodynamic Diameter	3 µm (MMAD) according to EN-13544-1
Respirable Fraction (< 5µm)	> 75%
Duty Cycle	30 min ON – 30 min OFF
Noise Level (EN 13544-1 Standard)	Approx. 50 dBA
Weight	0,46 kg
Size (Main Unit)	cm 14 x 3 x 9
Years of Warranty	5
Place of Manufacturing	Italy
Shipping carton	6

CE 0123