

MIKO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with compact size and carrying bag for easy transport of unit, accessories and medications. Recently upgraded for non-stop operation in case intensive treatments are required.

MIKO ref.
(with carrying bag)

RE 300600/03

MIKO BASIC ref.
(without carrying bag)

RE 300600/12



STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
 - Adult Masks with elastic strings
 - Child Masks with elastic strings
 - Mouthpiece
 - Nosepiece (non invasive)
 - Air Tube
 - Air Filter
 - Carrying Bag (MIKO model only)
- For spare consumables see page 27*



TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	170 VA
Max Pressure	2.2 bar 32 psi 220 kPa
Operating Pressure	0.8 bar 12 psi 80 kPa
Max Air Flow	12 l/min
Operating Air Flow	4,8 l/min
Neb-rate with 4 ml solution	0,35 ml/min
Mass Median Aerodynamic Diameter	2,44 µm (MMAD) according to EN-13544-1
Respirable Fraction (< 5µm)	> 80%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 56 dBA
Weight	1,45 kg
Size (Main Unit)	cm 14 x 17 x 11
Years of Warranty	3
Place of Manufacturing	Italy
Shipping carton	6

CE 0123